

High Performance Dioxin/Furan Immunoassay Kit

Materials List for DF1-12 Kit (for screening analysis of 12 samples)

- Kit insert (IN-DF1)
- Brief Summary of Immunoassay Protocol (PS-DF1)
- 20 Anti-Dioxin “**Antibody Coated Tubes**” (1 bag of 20 tubes with indicating desiccant)
- 1 bottle “**Sample Diluent**” (12 mL/bottle)
- 1 bottle “**Competitor-HRP Conjugate**” (12 mL/bottle)
- 1 bottle “**HRP Substrate Solution**” (12 mL/bottle)
- 1 bottle “**Stop Solution**” (12 mL/bottle)
- 1 vial of 0.5 mL of neat **Triton X-100**
- 1 vial of “TEG-Triton Keeper Stock” solution (Triton X-100 in tetraethylene glycol)

