

High Performance Dioxin/Furan Immunoassay Kit

Materials List for Kit DF1-60 (for screening analysis of 60 samples)

- Kit insert (IN-DF1)
- Brief Summary of Immunoassay Protocol (PS-DF1)
- 100 Anti-Dioxin “**Antibody Coated Tubes**” (5 bags of 20 tubes with indicating desiccants)
- 1 bottle “**Sample Diluent**” (55 mL/bottle)
- 1 bottle “**Competitor-HRP Conjugate**” (55 mL/bottle)
- 1 bottle “**HRP Substrate Solution**” (55 mL/bottle)
- 1 bottle “**Stop Solution**” (55 mL/bottle)
- 1 vial of 0.5 mL of neat **Triton X-100**
- 2 vials of “TEG-Triton Keeper Stock” solution (Triton X-100 in tetraethylene glycol)

