



# Immunoassay Sample Preparation Starter Kit SP2-ST

## Intended Purpose

for use with Sample Preparation Kits labeled SP2 and higher

contains non-disposable items for cleanup of 6 extracts using manually pressurized column method

## Materials List

- 7 glass reservoirs for carbon column wash
- 6 sets of stopper/stopcock for pressurizing columns
- 2 polypropylene syringes for pressurizing columns (10 mL & 20 mL)
- 2 Luer plugs to fit carbon columns for vacuum depressurization