



## Immunoassay Sample Preparation Kit

### Materials List for Kit SP3-60

**(for extraction and cleanup of 60 soil or sediment samples)**

- 60 wooden spatulas (5 bags of 12)
- 60 extraction vials (2 boxes of 30)
- 180 stainless steel mixing beads (5 bags of 36)
- 5 vials of bulk acid silica for extract pretreatment
- 60 individually sealed acid silica columns (5 boxes of 12)
- 60 activated carbon mini-columns (5 moisture proof pouches of 12, with indicating desiccant)